FOAM 1 & 2 CLASSIC

FOAMING CHARACTERISTICS OF LUBRICATING OILS TESTER



STANDARDS

ASTM D 892, D 6082, IP 146, ISO 6247 and related methods.

SC07E

This test method covers the determination of the foaming characteristics (foam tendency and stability) of lubricating oils at 24°C and 93.5°C.



*Photos are non-contractual - they may differ from actual products.



TECHNICAL FEATURES

- Test bath at 24°C and/or 93.5°C
- Fitted with cooling coil for test under the ambient temperature and better stabilization at 24°C
- 2 test position / bath
- Electronic temperature regulation with digital display
- Calibrated buit-in flow meter for air flow adjustment and control (2 / bath)
- Easy coupling for inlets/outlets (colored)
- FOAM 1: for sequences I and II
- FOAM 2: for sequences I, II and III

ACCESSORIES

- NDW Tech, diffuser washer with integrated pump (P/N: 941645)
- Wet test gas meter (P/N: 941644)
- Cooling cryostat (different model available depending on your application)
- Air compressor (P/N: 9426433)
- Drying tower ready for use (P/N: 9416401)
- Stainless steel spherical diffuser (different model available depending on your application)

SCOPE OF DELIVERY

941640 - FOAM 1

FOAM 1 is delivered with:

- 2 spherical diffusers (P/N: 24795)
- 2 test cylinders
- 2 flowmeters
- 1 ASTM 12C thermometer
- 1 bath at 24°C and 93.5°C

Site requirements:

- Power supply: 230V 50 Hz 10A
- Dimensions: (W) 340 x (D) 340 x (H) 770mm
- Weight: 20kg

941643 - FOAM 2

FOAM 2 is delivered with:

- 2 test baths at 24°C and 93.5°C
- 2 ASTM thermometer (12C), range from -20°C to +102°C (without certificate)
- 4 diffusers

Site requirements:

- Power supply: 230V 50 Hz 10A
- Dimensions: (W) 880 x (D) 340 x (H) 770mm
- Weight: 40kg



CONTACT: sales@normalab.com

Normalab FRANCE SAS 175, rue Claudie HAIGNERE 76190, Valliquerville

Tel.: +33 232.700.100 Fax: +33 232.704.732 DISTRIBUTED BY